

Z. klin. Chem. u. klin. Biochem.
10. Jg. 1972, S. 220—225

Useful Buffer and Gel Systems for Polyacrylamide Gel Electrophoresis

By H. R. MAURER and R. C. ALLEN

Max-Planck-Institut für Virusforschung, Tübingen, and Institut für gerichtliche Medizin der Universität Heidelberg

(Eingegangen am 10. November 1971)

Seven buffer and gel systems for polyacrylamide gel electrophoresis using discontinuous voltage and pH gradients (disc electrophoresis according to ORNSTEIN-DAVIS) and five systems using discontinuous voltage gradients at constant pH (according to ALLEN) have been compiled and characterized. Usefulness and applicability of the systems for clinical investigations are discussed and illustrated by several serum protein separations.

Sieben bewährte Puffer- und Gelsysteme für die Polyacrylamidgel-Elektrophorese mittels diskontinuierlicher Spannungs- und pH-Gradienten (Disk-Elektrophorese nach ORNSTEIN-DAVIS) und fünf Systeme für die Gel-Elektrophorese mittels diskontinuierlichem Spannungsgradienten bei konstantem pH-Wert (nach ALLEN) wurden zusammengestellt und charakterisiert. Die Brauchbarkeit der Systeme wurde diskutiert und an mehreren Beispielen von Serumprotein-Trennungen erläutert.

Although the high resolving power of polyacrylamide gel electrophoresis has proved valuable in clinical, biochemical and biological research at the molecular level (1—3), the usefulness and applicability of this method has been relatively little recognized or utilized in practical clinical chemistry for routine diagnosis and therapy control. There may be two main reasons for this.

First, there are, in some cases, many possible interpretations of the complicated protein patterns obtained with polyacrylamide gel electrophoresis. At the moment, a satisfactory interpretation of all the complex patterns is not possible. An example of these difficulties in interpretation is the overlapping of α_2 -globulins with β - and γ -globulins of normal human serum. The γ -globulins represent a broad complex which may be difficult to evaluate, since it extends as a protein background for α_2 - and β -globulins from the start to the α_2 -zone. Suitable reference systems, for the identification and characterization of the proteins under study may help to solve these method-inherent problems.

A second obstacle for the wide use of the method may be a technical one. Although there are various suitable apparatus systems commercially available, it appears that the techniques involved in the operation of most apparatus still require a certain degree of practical skill to obtain reproducible results. Moreover, as yet more correlative studies between laboratories are needed to provide satisfactory data on the standardization of sample preparation, separation and staining as well as on the standardization of gel formation. To this end, a further requirement is for buffer and gel systems which through adequate testing have been found to yield reproducible separations among different laboratories.

We have compiled several buffer and gel systems which have been successfully used in our and other laboratories. Although a sophisticated theory of multi-

phasic buffer systems has been formulated and a computer program based on this theory developed (4), we trust that the buffer and gel systems recorded in our tables will yield satisfactory separations in many cases, since they have been practically examined many times and found to work efficiently. These systems will be described and their applicability discussed. Examples illustrating the type of separations obtainable with various buffer and gel formulations will be presented. Problems related to the standardization of polyacrylamide gel electrophoresis for use in clinical chemistry will be dealt with in a subsequent paper (5).

Characterization of and Comments on the Gel and Buffer Systems for Polyacrylamide Gel Electrophoresis

Table 1 lists several gel and buffer systems based on the principles of the original disc electrophoresis as formulated by ORNSTEIN and DAVIS (6, 7). This technique utilizes both a discontinuous voltage and pH gradient to create the conditions of the KOHLRAUSCH regulating function for zone sharpening of sample components. While a continuous system employs the same buffers for electrode reservoirs and supporting medium (gel), a discontinuous system is composed of two buffers with a slower moving (trailing) ion in the electrode buffer and a faster moving (leading) ion in the gel and sample buffer. During electrophoresis the trailing ion displaces the position occupied by the leading ion, thereby sandwiching the sample components at the moving boundary. The principles of disc electrophoresis require pH differences between electrode, spacer and separation gel buffer to achieve the needed effective mobilities of the trailing ions. It is clear that the rather broad pH shifts during disc electrophoresis may, in some cases, give rise to concentration artefacts (irreversible protein interactions) and to the loss of bio-

Wächst Ihrem Labor die Arbeit über den Kopf?


Mit dem SILAB-System können Sie den gesamten Untersuchungsablauf praktisch voll automatisieren – von der Untersuchungsanforderung über die Probenanalyse bis zur Ergebnisdarstellung. Das SILAB-System ist nach dem Baukastenprinzip in drei Einzelsysteme gegliedert:

- Das SILAB-Identifikationssystem ordnet Untersuchungsgut und Meßdaten dem Patienten unverwechselbar zu.
- Das SILAB-Datenerfassungssystem sammelt, normiert und konzentriert die von den verschiedenen Laborgeräten anfallenden Daten.
- Das SILAB-Programmsystem, ein Software-Paket, bietet Programme für laufende Meßwertkontrolle, automatische Ergebnisberechnung und individuelle Ergebnisausgabe.

Und noch ein entscheidendes Argument für SILAB: Ob in Klinik oder Praxis – Sie können mit unserem System medizinische Laboratorien jeder Größenordnung voll- oder teilautomatisieren.


Für weitere Informationen wenden Sie sich bitte an unsere nächstgelegene Geschäftsstelle oder an:

SIEMENS AKTIENGESellschaft
Med Dv · 8520 Erlangen · Postfach 400



Technischen Fortschritt nützen mit dem SILAB-System von Siemens

Medizinischer DOLMETSCHER in Bildern



von Professor Dr. med.

Jack Adams-Ray

Verständigung leicht gemacht — im Ausland, mit Gastarbeitern und ausländischen Touristen

über alle Probleme eines Kranken mit Arzt und Schwester


durch anschauliche Symbole. Kurztexte deutsch, englisch, französisch, griechisch, italienisch, polnisch, russisch, serbo-kroatisch, spanisch, türkisch

**Einzelstück DM 5.— (iPr)
ab 10 Stück DM 4.—, ab 100 DM 3.—**



J.F. LEHMANNS VERLAG MÜNCHEN

Ärger mit den Schläuchen? Wir bei Heidolph wissen darum und bieten Ihnen deshalb die Pumpe SP-V an.



Nicht mehr starr der Abstand zwischen Anpreßflächen und Rollen, sondern anpassungsfähig an Schläuche mit Wandstärken zwischen 0,5 und 2,5 mm und lichten Weiten zwischen 2 und 8 mm. Und das ohne Werkzeug!

Stufenlos regelbar die Drehzahl. Abzulesen auf dem angebauten Drehzahlmesser.

Damit Sie dies alles als wesentliche Vorteile beurteilen können, senden wir Ihnen auf Anforderung gerne ausführliches Informationsmaterial.

Heidolph

Heidolph-Elektro KG · D 8420 Kelheim/Germany
Telefon (0 94 41) 75 33 · Telex 065415d

Gel systems for polyacrylamide gel electrophoresis using discontinuous voltage and pH gradients (disc electrophoresis) according to ORNSTEIN-DAVIS

No.	Designation	Gel system			Separation gel stock solutions			No.	Spacer gel stock solutions			Electrode buffer solution		Polarity	
		Range of application (examples)	Conc. at pH	Separation at pH	Components per 100 ml solution	pH	Mixing ratio (v/v)		Components per 100 ml solution	pH	Mixing ratio (v/v)	Components per 100 ml solution	pH	Top	Bottom
1	pH 8.9—7.0% (medium pore gel) (6)	Proteins of MW 10^4 — 10^6 ; optimum resolving range from $3 \cdot 10^4$ — $3 \cdot 10^5$ (e.g. serum proteins)	8.3	9.5	1 Tris 36.6 g to pH 8.3 (~48 ml) 1 N HCl 2 TEMED 0.23 ml 2 Acrylamide 28.0 g 3 Bis 0.735 g 3 Per 0.14 g	8.9	1 p. No. 1 2 p. No. 2 1 p. H ₂ O 4 p. No. 3	4	1 Tris 5.98 g to pH 6.9 (~48 ml) 1 N HCl 5 TEMED 0.46 ml 5 Acrylamide 10.0 g 6 Bis 2.5 g 6 Riboflavin 4.0 mg 7 Sucrose 40.0 g	6.9	1 p. No. 4 2 p. No. 5 1 p. No. 6 4 p. No. 7	Tris 6.0 g Glycine 28.8 g For use: 10% aq. soln.	8.3	—	+
2	pH 7.9—7.0% (medium pore gel) (7)	Proteins unstable above pH 8.5	7.3	8.5	1 Imidazole 11.7 g to pH 7.9 1 N HCl 8 TEMED 0.24 ml 1 N HCl	7.9	1 p. No. 8 2 p. No. 2 1 p. H ₂ O 4 p. No. 3	9	1 Imidazole 2.93 g to pH 5.9 1 N HCl 1 N HCl	5.9	1 p. No. 9 2 p. No. 5 1 p. No. 6 4 p. No. 7	Imidazole 0.465 g DL-asparagine 14.4 g For use: 10% aq. soln.	7.3	—	+
3	pH 7.5—7.0% (medium pore gel) (8)	Proteins (especially enzymes) showing optimum separation at pH 8.0 and stability above pH 8.0	7.0	8.0	1 Tris 6.85 g to pH 7.5 (~48 ml) 1 N HCl 10 TEMED 0.46 ml 1 N HCl	7.5	1 p. No. 10 2 p. No. 2 1 p. H ₂ O 4 p. No. 3	11	1 Tris 4.95 g to pH 5.5 1 N H ₃ PO ₄ 1 N HCl 1 N HCl	5.5	1 p. No. 11 2 p. No. 5 1 p. No. 6 4 p. No. 7	Diethylbarbituric acid 5.52 g Tris 1.0 g For use: 10% aq. soln.	7.0	—	+
4	pH 7.3—7.5% (medium pore gel) (9)	Basic proteins migrating back to cathode in gel system No. 1	8.3	6.6	1 N KOH 8.0 ml Glycine 19.0 g 12 TEMED 0.08 ml 13 Acrylamide 60.0 g 14 Bis 1.6 g 14 Per 0.56 g	7.3	6 p. No. 12 1 p. No. 13 1 p. No. 14	15	1 N KOH 48.0 ml Glycine 4.8 g 1 N HCl 1 N HCl	10.3	1 p. No. 15 2 p. No. 5 4 p. No. 7 1 p. No. 14	Glycine 13.7 g 2,6-dimethylpyridine (2,6-lutidine) 38.2 ml For use: 10% aq. soln.	8.3	+	—
5	pH 7.2—7.0% (10)	S _α , S _β and γ-globulins	8.6	8.6	1 Tris 10.65 g to pH 7.2 1 N HCl 1 N HCl	7.2	1 p. No. 16 2 p. No. 2 1 p. H ₂ O 4 p. No. 3	17	1 Tris 6.0 g to pH 7.3 1 N HCl 1 N HCl	7.3	1 p. No. 17 2 p. No. 5 1 p. No. 6 4 p. No. 7	Tris 10.0 g Glycine 28.8 g For use: 10% aq. soln.	8.5	—	+
6	pH 4.3—15% (small pore gel) (11)	Basic proteins (e.g. histones) of MW about $2 \cdot 10^4$	5.0	3.8	1 N KOH 48.0 ml AcOH 17.2 ml 18 TEMED 4.0 ml 19 Acrylamide 60.0 g 20 Bis 0.4 g 20 Per 0.28 g	4.3	1 p. No. 18 2 p. No. 19 1 p. H ₂ O 4 p. No. 20	21	1 N KOH 48.0 ml AcOH 2.87 ml 1 N HCl 1 N HCl	6.7	1 p. No. 21 2 p. No. 5 1 p. No. 6 4 p. H ₂ O	β-Alanine 31.2 g AcOH 8.0 ml For use: 10% aq. soln.	4.5	+	—
7	pH 2.9—7.0% (medium pore gel) (7)	Highly basic proteins of MW $\sim 2 \cdot 10^4$	4.0	2.3	1 N KOH 12.0 ml AcOH 53.25 ml 22 TEMED 1.15 ml 23 Per 2.8 g	2.9	4 p. No. 22 2 p. No. 2 2 p. No. 23	24	1 N KOH 48.0 ml AcOH 2.95 ml 1 N HCl 1 N HCl	5.9	1 p. No. 24 2 p. No. 5 1 p. No. 6 4 p. H ₂ O	Glycine 28.1 g AcOH 3.06 ml For use: 10% aq. soln.	4.0	+	—

The designation "pH 8.9—7.0%-gel" means that the separation gel of the corresponding gel system contains a buffer of pH 8.9 and 7.0% acrylamide (w/v). "Concentration" and "separation" at pH 8.3 and 9.5, respectively, indicate at which pH values the substances are concentrated (stacked) and separated during electrophoresis. All solutions should contain double-distilled water as solvent. The polarity is applicable with the assumption that the sample components migrate vertically downward. Distilled water can be used in place of the saccharose (sucrose) solution (No. 7). The molecular weight values are approximate. All ammonium persulfate solutions should be freshly prepared each day and should be added to the other solutions just before pouring the gel. All solutions should be at room temperature (20°) before mixing.

Abbreviations:

AcOH Glacial acetic acid
Bis N,N'-methylene-bis-acrylamide
EDMA Ethylenediacrylate
MW Molecular weight
p. Volume parts

Ammonium persulfate (NH₄)₂S₂O₈
Sodium bisulfite, NaHSO₃
N,N,N',N'-tetramethylethylenediamine
Tris(hydroxymethyl)aminomethane as base separation

Per SO₃
TEMED Tris
sep. sep.

Tab. 2
Gel systems for polyacrylamide gel electrophoresis using discontinuous voltage gradients (ALLEN-System)

No.	Designation	Range of application (examples)	Gel stock solutions		pH at 25°	Mixing ratios (v/v) of stock solutions for step gradient gels				Electrode buffer solution	pH at 25°
			No.	Components per 100 ml solution		1st. sep. gel	2nd sep. gel	3rd sep. gel	cap and well forming gel		
1	pH 9-plasma-protein-system	Plasma and serum proteins, esterases, tissue proteins	1 (1.5 M) 2 3 3a 3b 4 (0.3 M)	1 N H ₂ SO ₄ Tris TEMED Acrylamide Bis Per Per Solin. No. 1 TEMED	31.0 ml 18.15 g 0.24 ml 32.0 g 0.75 g 0.21 g 0.105 g 70 mg 20.0 ml 0.19 ml	8% gel: 10.0 p. No. 1 10.0 p. No. 2 13.5 p. H ₂ O 6.5 p. No. 3 Sample solution (No. 5, 0.075 M): Solin. No. 1 (less TEMED) Sucrose Bromphenolblue, 0.1% aq. H ₂ O to	6% gel: 12 p. No. 1 9 p. No. 2 19 p. H ₂ O 8 p. No. 3	3.5% gel: 8.0 p. No. 1 3.5 p. No. 2 14.5 p. H ₂ O 6.0 p. No. 3	8% gel: 1 p. No. 4 1 p. No. 2 2 p. No. 3 pH 9.0	Tris (0.065 M) Boric acid Methiolate (Thiomersol) 0.1 g	9.0
2	pH 9-lipo-protein-system	Plasma lipo-proteins	6 (0.75 M) 7 (0.15 M)	1 M Citric acid Tris TEMED Brij 35, 10% aq. Solin. No. 6 TEMED Brij 35, 10% aq.	7.0 ml 9.075 g 0.24 ml 0.2 ml 20.0 ml 0.19 ml 0.2 ml	8% gel: 1 p. No. 6 1 p. No. 2 2 p. No. 3a Sample solution (No. 8, 0.0375 M): Solin. No. 6 (less TEMED) Sucrose Bromphenolblue, 0.1% aq. H ₂ O to	5% gel: 16 p. No. 6 10 p. No. 2 6 p. H ₂ O 32 p. No. 3a	3.5% gel: 8.0 p. No. 6 3.5 p. No. 2 14.5 p. H ₂ O 6.0 p. No. 3	8% gel: 1 p. No. 7 1 p. No. 2 2 p. No. 3 pH 9.0	Tris (0.065 M) Boric acid Methiolate (Thiomersol) 0.1 g	9.0
3	pH 8.5-enzyme-system	Spleen extract proteins, esterases, LDH isoenzymes	9 (1.5 M) 10 (0.3 M)	1 N HCl Tris TEMED Solin. No. 9 TEMED	58.0 ml 18.15 g 0.24 ml 20.0 ml 0.19 ml	8% gel: 1 p. No. 9 1 p. No. 2 2 p. No. 3a Sample solution (No. 11, 0.075 M): Solin. No. 9 (less TEMED) Sucrose Bromphenolblue, 0.1% aq. H ₂ O to	6% gel: 12 p. No. 9 9 p. No. 2 3 p. H ₂ O 24 p. No. 3a	3.5% gel: 8.0 p. No. 9 3.5 p. No. 2 14.5 p. H ₂ O 6.0 p. No. 3	8% gel: 1 p. No. 10 1 p. No. 2 2 p. No. 3 pH 8.5	Tris Glycine 6.0 g 28.75 g	8.5
4	pH 9-nucleic acid-system	Nucleic acids, 16-4 S RNA	12 13	Acrylamide Bis Solution No. 1 TEMED	48.0 g 1.68 g 50.0 ml 0.12 ml	12% gel: 1 p. No. 12 1 p. No. 1 2 p. No. 3b Sample solution (No. 14, 0.187 M): Solin. No. 1 (less TEMED) Ethylene glycol Bromphenolblue, 0.1% aq.	7.5% gel: 3 p. H ₂ O 5 p. No. 12 8 p. No. 1 16 p. No. 3b	4.5% gel: 7.5 p. H ₂ O 4.5 p. No. 12 12.0 p. No. 1 24.0 p. No. 3b	8% gel: 3.4 p. H ₂ O 6.6 p. No. 12 4.6 p. No. 13 10.0 p. No. 13 20.0 p. No. 3 pH 9	As for System No. 1	9.0
5	pH 4-basic protein-system	Basic proteins, histones, ribosomal proteins	15 (0.48 M) 16	1 N KOH Glacial acetic acid TEMED 1 N KOH Glacial acetic acid TEMED	48.0 ml 27.0 ml 4.0 ml 12.0 ml 14.0 ml 4.0 ml	12% gel: 1 p. No. 12 1 p. No. 15 2 p. No. 3 Sample solution (No. 17, 0.03 M of K ⁺): Sucrose 0.12 M Potassium acetate pH 4 H ₂ O Final volume	8% gel: 1 p. No. 2 1 p. No. 15 2 p. No. 3	4.5% gel: 8.0 p. No. 15 4.5 p. No. 2 3.5 p. H ₂ O 16.0 p. No. 3a	8% gel: 1 p. No. 2 1 p. No. 16 2 p. No. 3 pH 4	β-Alanine Glacial acetic acid Dilute 1:20 for use 58 g 68 ml	4.0

The samples are dissolved in or diluted with the corresponding sample solution. The polarity is: cathode ⊖ top, anode ⊕ bottom, except for system No. 5 (reverse). Staining see Appendix. Abbreviations: see note to Table 1

logical activity. Therefore, attention should always be given to pH influences to which the sample components are exposed in every disc electrophoretic run. For example, in the standard alkaline gel system No. 1 of Table 1, the sample components are dissolved in spacer gel buffer of pH 6.7, are concentrated in the spacer gel at pH 8.3 and are separated in the separation gel at pH 9.5.

Moreover, the ion species used for the moving boundary may influence the activity of enzymes. E. g., the alkaline phosphatase of serum which shows several bands in polyacrylamide gel electrophoresis is inhibited by glycinate, the trailing ion of gel system No. 1. Inactivation of enzymes may also be caused by photopolymerization of the sample gel or by interaction with other catalysts of gel polymerization. The pertinent aspects of polymerization artifacts are discussed elsewhere (1).

Gel system No. 1 is widely used for the separation of proteins (e. g. of serum) of molecular weights between 10^4 and 10^6 which are stable at the pHs indicated (7). System No. 2 (8) and No. 3 (9) are suited for proteins which are unstable above pH 8.5 and 8.0, respectively. System No. 4 (10) catches basic proteins which migrate retrograde (to cathode) in gel system No. 1. The latter gel system normally separates anionic substances with free mobilities in the range of -0.6 to -7.5 units. However, proteins with intermediate mobilities (-0.75 to 1.0), such as the S_{α} -, S_{β} - and γ -globulins, may be better resolved by gel system No. 5 (11). Gel systems No. 6 (12) and 7 (8) are suited for basic proteins, such as histones, proteohormones and ribosomal proteins. To create a steeper voltage gradient, resulting in improved resolution, glycine may be replaced by valine in the anionic (upper) electrode buffer of gel system No. 7. Valine has a similar buffering capacity to glycine, but a lower mobility at pH 4.0 (13).

The gel and buffer systems of Table 2, developed by ALLEN (14–16), separate the sample components on the basis of discontinuous voltage gradients at a uniform pH throughout. Essentially, they combine two methods for sample concentration, i. e., the principle of conductivity shift described by HJERTÉN (17) and the moving boundary technique first described by POULIK (18). A conductivity shift is readily generated by lowering the ionic strength of the sample buffer with respect to that of the separating medium. This concentrating method avoids possible harmful effects arising from pH discontinuities, but has a limited concentrating capacity. The moving boundary technique, operating under the conditions described by the KOHLRAUSCH regulating function, sharpens the zones of already separated components, thus improving resolution considerably. The technique also uses discontinuous buffers with leading and trailing ions, which, however, in contrast to the disc electrophoresis systems, have a uniform pH throughout the total system, since the mobility of the trailing ion is constant in this system.

Gel system No. 1 of Table 2, a sulfate-borate discontinuous buffer system, has proved valuable for the separation of many serum proteins, enzymes and tissue proteins which are stable at pH 9.0. In order to utilize the potential of the molecular sieving effect of the gel, a step gradient gel is recommended: First a 8% gel (with respect to acrylamide), then a 6% and a 3.5% gel are polymerized on top of each other. For Gc type proteins a 12 — 8 — 6 — 3.5% step gradient gel is advised. Following sample application into the pockets of the well gel, an additional cap gel is cast on top of the samples to physically separate the samples from the trailing ions in the upper electrode reservoir. Thus, the moving boundary formed by the leading-trailing ions is maintained above the sample components during initial separation by charge and size in the first separating gel. The ionic strength of the cap gel is less than that of the separating gel to increase the velocity of the moving boundary. If the persulfate used as catalyst for the gel formation interferes with the enzyme or other activity of the biological material to be separated, the well and cap gel may be photopolymerized rather than chemically polymerized.

Addition of thioglycolate or mercaptoethanol (5 mM) may help to reduce the oxygen-rich state of the gel due to persulfate decomposition (19, 20).

The citrate-borate system No. 2 is suited for serum proteins, particularly lipoproteins, which tend to aggregate. Citrate and borate are chelating agents which might prevent extensive aggregation due to ion interactions. Addition of a non-ionic detergent, such as Brij 35, supports the dissociation of clustering lipoproteins. The tris-glycine system No. 3 is recommended as a substitute for the borate systems, if this ion should complex with sugar and carbohydrate-containing components leading to bent protein zones (21). Gel system No. 4 is a modified gel system No. 1 for optimum separation of nucleic acids (4–16 S RNA). Finally, gel system No. 5 was developed on the basis of gel system No. 6 of Table 1 for the fractionation of basic proteins at pH 4.

It should be added that reducing the ionic strength of the separating gels of systems No. 1 and 3 to one half has recently been found to give equal or better resolution, utilizing lower voltage and lower current densities with reduced Joule heating.

Figure 1 shows the type of pattern of human serum proteins obtainable with the standard ORNSTEIN-DAVIS system No. 1 (Tab. 1), a discontinuous glycine-chloride buffer system, while Figure 2 demonstrates human serum patterns produced with ALLEN system No. 1, a pH-continuous, sulfate-borate-discontinuous buffer system using a 8 — 6 — 3.5% step concentration gradient. A comparison reveals that the latter system yields more and sharper bands, particularly in the pre- and postalbumin region. This difference is more pronounced, hence the resolution is improved, if the ionic strength of ALLEN's system is reduced to one half and the step gradient changed to include a 12% gel (Fig. 3).




Fig. 1
Human serum proteins in standard ORNSTEIN-DAVIS gel system No. 1 (Tab. 1) in a flat slab gel (3 mm thick). Sample quantities: 200 μ g protein on each. Constant current of 60 mA. Coomassie Blue R 250 staining

⊕

⊖




Fig. 2
Human serum proteins in ALLEN's gel system No. 1 (Tab. 2) in a flat slab gel (3 mm thick) using a 8—6—3.5% step concentration gradient. Other conditions see Figure 1

⊕

The improved separation obtained by ALLEN's system is mainly due to the step (discontinuous) polyacrylamide gradient which sharpens boundaries as they pass gel concentration discontinuities. Gradient gel electrophoresis considerably increases the versatility of polyacrylamide gel electrophoresis and deserves, therefore, particular attention if high resolution of distinct protein bands is desired. The method includes several parameters, and hence it offers a wide choice of conditions: Different acrylamide and/or cross-linker concentration gradients can be formulated; moreover, the form (linear, stepped, concave, convex) and the slope of the gradient may be modified. Methods for the preparation of gradient gels and problems related to

these methods are described and discussed elsewhere (1).

It follows from theoretical considerations and practical observations that there hardly exists a single gel system which would be optimal for the resolution of all proteins of a complex biological fluid. It is rather advisable, in developing a particular gel system, to limit the goal in each case to optimum separation of the proteins of interest. Two equations developed by RODBARD and CHRAMBACH (2, 22, 23) may serve as guidelines for the optimization of separation and resolution. These require knowledge of the retardation coefficient and the free mobility to calculate the gel concentration for maximal separation and for optimal resolution. More-




Fig. 3

Human serum proteins in ALLEN's gel system No. 1 (Tab. 2) using half ionic strength, a 12—8—6—4.5% step concentration gradient and pulsed constant power (280 V, 1.0 MFD, pulse rate increased from 75—300 PPS over first 20 min. in 3 steps of 75 pulses each at 5 min. intervals). Alternating normal serum and myeloid leukemia serum proteins (from left to right) separated. Other conditions see Figure 1

over, JOVIN and coworkers (3, 4) developed a theory of multiphasic buffer systems and a computer program based on this theory permitting operation at any pH, at 0° or 25°, thus introducing a new dimension of

versatility into fractionation. For example, electrophoresis at several pH values allows charge fractionation, the construction of pH activity and stability profiles as well as titration curves. The usefulness and applicability of these systems in clinical chemistry have still to be examined.

Appendix

New Staining Procedures

Proteins with Coomassie Brilliant Blue R 250

The proteins are first fixed by incubating the gel in 12.5% trichloroacetic acid at 65° for 30 min. The gel is rinsed with tap water and stained with a mixture of Coomassie Brilliant Blue R 250 (0.2% water/absolute ethanol/glacial acetic acid, 45:45:10, V/V) at 65° for 30 min. The gel should be turned over after 10 min. in the stain. This procedure is also suited for isoelectric focused gels.

Lipoproteins with Lipid Crimson

3 g of Lipid Crimson (see (1)) are dissolved in 600 ml of 95% ethanol, 400 ml water are added and the solution allowed to age one week. Following filtration, the gel is stained 2 hrs at 65° and left overnight at room temperature in the same stain. Destaining is performed with 10% acetic acid, with 3 to 4 changes in a 24 hr period.

RNA with Stains-All (25)

Formamide is adjusted to pH 7.3—7.4 with concentrated HCl. 0.1% of the dye is dissolved and the stock solution stored protected from light at 4°. The working solution is freshly prepared by mixing 5 volume parts of stock solution, 45 parts of formamide (pH 7.3) and 50 parts of water. The gel is stained 1 hr at room temperature in a dark place. Destaining is with water. The stain rapidly fades when exposed to strong light.

References

- MAURER, H. R., Disc Electrophoresis and Related Techniques of Polyacrylamide Gel Electrophoresis, Walter de Gruyter, Berlin-New York, 224 pages, 2nd revised and expanded edition 1971. First edition: Disk-Elektrophorese, Theorie und Praxis der diskontinuierlichen Polyacrylamidgel-Elektrophorese, Walter de Gruyter, Berlin 1968. — 2. GORDON, A. H., in: Laboratory Techniques in Biochem. and Mol. Biol. Vol. 1, part 1, Work, T. S., and Work, E., editors, North-Holland Publ. Co. p. 1—145, Amsterdam-London (1969). — 3. CHRAMBACH, A. and D. RODBARD, Science 172, 440 (1971). — 4. JOVIN, T. M., M. L. DANTE and A. CHRAMBACH, Multiphasic Buffer Systems Output (Nat. Techn. Inform. Serv. Springfield, Va., USA) PB 196085 to 196091 (1970). — 5. MAURER, H. R. and R. C. ALLEN, in preparation. — 6. ORNSTEIN, L., Annals N. Y. Acad. Sci. 121, 321 (1964). — 7. DAVIS, B. J., Annals N. Y. Acad. Sci. 121, 404 (1964). — 8. Canalco Chemical Formulations for Disc Electrophoresis (1965). — 9. WILLIAMS, D. E. and R. A. REISFELD, Annals N. Y. Acad. Sci. 121, 373 (1964). — 10. TABER, H. W. and F. SHERMAN, Annals N. Y. Acad. Sci. 121, 600 (1964). — 11. JOLLEY, W. B. and H. W. ALLEN, Nature 208, 390 (1965). — 12. REISFELD, R. A., U. J. LEWIS and D. E. WILLIAMS, Nature 195, 281 (1962). — 13. SHEPHERD, G. R. and L. R. GURLEY, Analytic. Biochem. 14, 356 (1966). — 14. ORTEC Instruction Manual Model 4200 Electrophoresis System (1969). — 15. ORTEC A 32 Techniques for High Resolution Electrophoresis (1970). — 16. ORTEC Instruction Manual Model 4217 Electrophoresis System (1971). — 17. HJERTÉN, S., S. JERSTEDT and A. TISELIUS, Analytic. Biochem. 11, 219 (1965). — 18. POULIK, M. D., Nature 180, 1477 (1957). — 19. BREWER, J. M., Science 156, 256 (1967). — 20. KING, E. E., J. Chromatog. 53, 559 (1970). — 21. LERCH, B. and H. STEGMANN, Analytic. Biochem. 29, 76 (1969). — 22. RODBARD, D. and A. CHRAMBACH, Proc. Nat. Acad. Sci., USA 65, 970 (1970). — 23. RODBARD, D., G. KAPADIA and A. CHRAMBACH, Analytic. Biochem. 40, 135 (1971). — 24. RODBARD, D. and A. CHRAMBACH, Analytic. Biochem. 40, 95 (1971). — 25. DAHLBERG, A. E., C. W. DINGMAN and A. C. PEACOCK, J. Mol. Biol. 41, 139 (1969).

Dr. Hans Rainer Maurer
Max-Planck-Inst. f. Virusforschung
74 Tübingen/Germany
Spemannstr. 35