

Z. Klin. Chem. Klin. Biochem.  
11. Jg. 1973, S. 184–188

## A Simple, Kinetic Procedure for the Optimized Determination of the Activity of Serum Glutamate Dehydrogenase (EC 1.4.1.2.) with an Enzyme Analyzer

By R. HAECKEL

*Institut für Klinische Chemie (Geschäftsführender Direktor: Prof. Dr. Dr. J. Büttner)  
Medizinische Hochschule Hannover*

(Eingegangen am 16. November/29. Dezember 1972)

The reliability of a mechanized procedure for the determination of the glutamate dehydrogenase activity in human sera was investigated with an Eppendorf analyzer 5011. In spite of a 1:3.5 ratio of sample to assay volume, carry-over effects can be neglected between sera containing less than 60 U/l. For routine purposes it is suggested that 0.3 U/l be subtracted from each activity found in human sera and that individual blank values are not measured.

Die Zuverlässigkeit einer mechanisierten Bestimmung der Glutamatdehydrogenase-Aktivität in menschlichem Serum wurde mit einem Eppendorf-Automaten 5011 untersucht. Obwohl das Verhältnis von Probe- zu Ansatzvolumen 1:3,5 beträgt, können Verschleppungseffekte zwischen Seren, die weniger als 60 U/l enthalten, vernachlässigt werden. Für Routineuntersuchungen wird vorgeschlagen, auf die Bestimmung eines individuellen Probenleerwertes zu verzichten und 0,3 U/l als konstanten Betrag von der im menschlichen Serum gefundenen Aktivität abzuziehen.

The activity of serum glutamate dehydrogenase (EC 1.4.1.3) has become a valid biochemical parameter for the differential diagnosis of hepatic diseases (1–3). For this purpose, the enzyme is determined routinely on a large scale at our hospital. Therefore, a mechanized procedure was evaluated for its reliability.

The activity of this enzyme was measured under "optimized" conditions according to SCHMIDT (4) and the recommendations of the German Society for Clinical Chemistry (5). This procedure includes the determination of a sample blank value. For purposes of simplification we have investigated whether this blank could be subtracted from the result as a fixed value.

### Materials and Methods

Triethanolamine, NADH (Trisodium salt) and ADP were purchased from C. F. Boehringer (D-68 Mannheim), control sera from Merz and Dade (D-8 Munich), bovine albumin (purest) from Behringwerke AG (D-355 Marburg) and all other chemicals p. a. from E. Merck AG (D-61 Darmstadt).

#### Solutions

1. Triethanolamine buffer: 90 mmol/l  
16.7 triethanolamine (Boehringer No. 15325), 1.675 g disodium-ethylenedinitriloacetate (Merck No. 8418) + bidist. H<sub>2</sub>O to 1000 ml. The pH value should be 8.0, otherwise it must be adjusted. This solution can be stored at 4°C for some months.
2. 2-Oxoglutarate: 252 mmol/l  
3.682 g 2-oxoglutarate (Merck No. 5194) + bidist. H<sub>2</sub>O to 100 ml. This solution can be stored 4 weeks at 4°C.
3. Ammonium acetate: 900 mmol/l  
6.937 g ammonium acetate (Merck No. 1116) + bidist. H<sub>2</sub>O to 100 ml. This solution is stable for several months at 4°C.

4. Nicotinamide-adenine-dinucleotide, reduced: 14.4 mmol/l  
51.1 mg NADH (Boehringer No. 15142) + bidist. H<sub>2</sub>O to 5 ml. This solution can be stored 24 h at 4°C.

5. Adenosine diphosphate: 72 mmol/l  
339.3 mg ADP (Boehringer No. 15691) + bidist. H<sub>2</sub>O to 10 ml (Stable for 1 week at 4°C).

6. Reaction mixture:

Solution 1	40 ml
Solution 3	8 ml
Solution 4	1 ml
Solution 5	1 ml

The reaction mixture should be prepared immediately before use; it can be stored for 6 h at 4°C.

#### Assay

Reaction mixture	500 $\mu$ l
Sample	200 $\mu$ l

Mix and incubate 15 min at 25°C

Solution 2	20 $\mu$ l
------------	------------

Mix and record the reaction at 334 nm

The end concentrations of the assay are the same as recommended by the German Society for Clinical Chemistry (5). Samples with an activity of glutamate dehydrogenase (GIDH) above 30 U/l are diluted 1 + 4 with Qualtrol (see Tab. 4).

Eppendorf micropipettes and an Eppendorf photometer (No. 1101 M + 2705)<sup>1)</sup> with an analogue printer (Eppendorf No. 4414) were used for the manual method. The mechanized procedure was performed with an enzyme analyzer 5011 from Eppendorf Gerätebau GmbH (D-2 Hamburg) following the instructions of the manufacturer (Programming plug: HBDH/LDH; filter: 334 nm; transformation: 0–0.5; chart speed: 1 cm/min). The glutamate dehydrogenase activity was calculated from the angle which is formed between the direction of the moving chart paper

<sup>1)</sup> Automatic cuvette-positioning attachment.




Fig. 1

The sequence of samples introduced into the Eppendorf enzyme automat 5011 for the determination of the glutamate dehydrogenase activity. The intersegment B can be repeated as often as required

and a line drawn through the 7 points printed out by the analyzer:  $\text{tg } \alpha \cdot 15 = \text{U/l}$ . The angle  $\alpha$  was measured with an Eppendorf scale No. 2780. A table can be prepared containing the calculated activity values ( $\text{tg } \alpha \cdot 15 - 0.3$ ) and the corresponding angle.

The sequence of samples introduced into the analyzer is shown in Figure 1. Drift effects can be neglected and need not be monitored because of the nature of the kinetic procedure. No drift of blank or standard (6) was observed in batches containing up to 60 samples.

The detection limit was derived from the interbatch variation ( $n = 15$ ) of water samples Nos 5, 6, 9 and 10 (Fig. 1) as the threefold standard deviation according to KAISER (7). The mean was calculated from all 4 values.

The interaction from low to high concentrations ( $Q_1$ ) and from high to low concentrations ( $Q_2$ ) was determined as recently described (8) and expressed as the percent interaction coefficient ( $Q = 100 \cdot q$ ; for the explanation of  $q$  see Fig. 1). A sequence of 5 identical samples (Enzatzol, lot No. 233A) with a relatively low activity of glutamate dehydrogenase (6 U/l) was followed by 5 identical sera with a high activity (110 U/l) and finally by 5 samples with the low activity. Sera with a high glutamate dehydrogenase activity were prepared by mixing 5 ml Enzatzol (Lot No. ET 233A) with 50  $\mu\text{l}$  diluted glutamate dehydrogenase suspension (2.5 ml Enzatzol + 10  $\mu\text{l}$  glutamate dehydrogenase No. 15140 from Boehringer-Mannheim). Monitrol II solutions containing various glutamate dehydrogenase activities were obtained similarly.

The temperature at  $25.0 \pm 0.1^\circ\text{C}$  was checked in the cuvettes with a tasto therm instrument (Braun, D-6 Frankfurt) calibrated against a thermometer.

For the determination of the precision 3 identical serum samples ( $x_a, x_b, x_c$ ) are analyzed daily (Fig. 1): the precision within-series ( $\text{CVs} = \frac{100 \cdot \bar{s}}{s}$ ) is calculated from the  $x_b$  — and  $x_c$  values

$$\left( \bar{x} = \frac{\sum_{i=1}^n (x_{bi} + x_{ci})}{2m} ; s = \sqrt{\frac{\sum_{i=1}^n (x_{bi} - x_{ci})^2}{2m}} \right),$$

the precision from-day-to-day from the  $x_b$  — (or the  $x_c$ ) values

$$\left( \bar{x} = \frac{\sum_{i=1}^n x_{bi}}{n} ; s = \sqrt{\frac{\sum_{i=1}^n (x_{bi} - \bar{x})^2}{n-1}} \right).$$

## Results and Discussion

The procedure described above for the mechanized determination of the glutamate dehydrogenase activity in human serum samples is precise enough for clinical diagnostic purposes (Tab. 1). Under these experimental conditions an activity of up to 30 U/l can be measured and up to 60 U/l if the recorder paper runs at 2 cm per minute. Samples with higher activities must be diluted 1 + 4. The detection limit was 0.3 U/l ( $\bar{x} = 0.26$ ;  $n = 4$ ).

The activity of the glutamate dehydrogenase was determined in serum samples from 97 patients by the automated and by a manual procedure. The data obtained with both methods show a good correlation (Fig. 2).

Since a relatively large sample volume is chosen, interaction effects have to be considered. The percent carry over from high to low activities ( $Q_2$ ) was  $1.5 \pm 0.4$

Tab. 1

The precision of the mechanized determination of the glutamate dehydrogenase (GIDH) activity. Monitrol II (3 ml) was mixed with 50  $\mu\text{l}$  enzyme suspension containing various glutamate dehydrogenase (GIDH) activities as described under methods

	$\bar{x}$	Precision within-series			$\bar{x}$	Precision from-day-to-day		
		(n)	s	CVs		(n)	s	CVd
Monitrol II (+ GIDH)	3.1	(20)	0.30	9.7	5.2	(19)	0.43	8.3
	6.5	(20)	0.30	4.6				
	11.4	(20)	0.40	3.5				
	14.1	(20)	0.39	2.8				
	26.7	(20)	0.51	1.9				
	32.4	(20)	0.49	1.5				
Enzatzol	5.2	(19)	0.17	3.3	5.2	(15)	0.44	8.4
	4.4	(19)	0.38	8.6*				
Pooled human serum	7.5	(20)	0.26	3.4				
	7.8	(20)	0.48	6.2*				

\* 100  $\mu\text{l}$  sample + 500  $\mu\text{l}$  reaction mixture; end concentration of the assay according to l. c. (4, 5)




Fig. 2

Comparison of glutamate dehydrogenase activities in various serum samples determined with automated and a manual procedure:  $y = 0.9666x + 0.0444$ ; coefficient of correlation  $r = 0.9942$  ( $n = 97$ )




Fig. 4

The distribution of serum blank values in samples from 1626 patients. Bidist.  $H_2O$  was substituted for 2-oxoglutarate




Fig. 3

The "interaction-safe range" (8) of the mechanized procedure for the determination of the glutamate dehydrogenase activity. The concentrations forming this range are calculated according to l. c. (8):  $C_{1,1}$  for the interaction from low to high concentrations and  $C_{2,1}$  for the interaction from high to low concentrations. The standard deviations were taken from Table 1. It is postulated that the error due to carry over effects should be lower than the threefold (resp. twice) value of the within-series precision (standard deviation)

determined on 10 different days, and could not be detected from low to high activities ( $Q_1$ ). Using a procedure recently developed (8), a so called "interaction-safe range" was calculated (Fig. 3) which demonstrated that carry over effects can be neglected between samples containing less than 60 U/l. If a sample must be diluted because of its high activity, the analysis of the following sample must also be repeated.

In the test procedure proposed by SCHMIDT (4) an individual blank value has to be determined in which the substrate (2-oxoglutarate) is replaced by distilled water. The distribution of 1626 blank values from various patients shown in Figure 4 is not Gaussian. The median is at 0.3 U/l. 93% of all values are not higher than 0.5 and 97% not higher than 0.8.

For routine purposes we suggest that 0.3 U/l be subtracted from each activity found in human sera and that the determination of individual blank values be omitted. However, it must be determined whether the discrimination at the borderline between normal and patho-

Tab. 2

Activity values (U/l) which are estimated falsely above the normal range by the proposed procedure among 1616 unselected human sera

Patient	A uncorrected activity	B as A, minus 0.3	C as A, minus individual blank value
1	4.0	3.7	1.9
2	4.0	3.7	2.9
3	3.5	3.2	1.9
4	3.7	3.4	2.1

Tab. 3

The influence of the sample volume on the glutamate dehydrogenase activity of various human sera. The volume of the reaction mixture (500  $\mu$ l) and the end concentrations of the assay (according to l. c. 5) were kept constant. A = assay value, B = sample blank value (bidist. H<sub>2</sub>O was substituted for 2-oxoglutarate), C = A-B. All figures are mean values (U/l) from 2 determinations

Sample volume	50 $\mu$ l			100 $\mu$ l			200 $\mu$ l			300 $\mu$ l		
No.	A	B	C	A	B	C	A	B	C	A	B	C
1	8.4	2.5	5.9	7.2	0.9	6.3	6.1	0.8	5.3	6.7	0.2	6.5
2	4.1	1.3	2.8	4.5	0.7	3.8	3.4	0.8	2.6	3.4	0.2	3.2
3	18.3	1.0	17.3	16.2	0.9	15.3	12.2	0.3	11.9	13.3	0.3	13.0
4	4.1	0.8	3.3	3.0	0.5	2.5	2.6	0.5	2.1	1.9	0.2	1.7
5	7.6	1.3	7.3	5.3	0.9	4.4	4.0	0.3	3.7	3.6	0.2	3.4
6	45.8	1.6	44.2	39.7	1.4	38.3	29.5	0.8	28.7	23.4	0.4	23.0
7	26.9	1.3	25.3	23.7	0.9	22.8	18.2	0.4	17.8	18.0	0.3	17.7
8	23.1	1.3	21.8	19.8	0.9	18.9	16.1	0.5	15.6	15.4	0.2	15.2
9	29.1	1.3	27.8	27.3	0.9	26.4	21.9	0.8	21.1	21.7	0.5	21.2
$\bar{x}$	18.6	1.3	17.3	16.3	0.9	15.4	12.7	0.6	12.1	11.9	0.3	11.6

logical is not altered by the procedure just suggested. This borderline is at 3 U/l according to SCHMIDT (4).

In 4 different sera out of 1626 the activity of glutamate dehydrogenase, calculated by subtracting 0.3, was above 3 U/l, whereas the results were below 3 U/l when the correct individual blank values were used (Tab. 2). Thus, the estimated activity of glutamate dehydrogenase was falsely above the normal range in 0.25% of all human sera checked during this investigation. In these cases, however, the error was too low to seriously mislead the physician. For scientific studies the individual blank values should further be considered.

Since a large serum volume is used for the glutamate dehydrogenase assay, the influence of sample volume on the enzyme activity was investigated. It was found with various sera that the glutamate dehydrogenase activity declines with increasing the sample volume (Tab. 3). This indicated that the proposed assay conditions are not optimal for serum samples. SCHMIDT (4) uses a 1:3.5 ratio between the volumes of the sample and the reaction mixture. If this ratio is changed to 1:6 and the end concentrations of the assay kept identical according to l. c. (5) the variation within series batch increased significantly (Tab. 1). Because of this result we suggest that the recommendation of SCHMIDT (4) be followed. This appears to be a reasonable com-

promise for routine analytical purposes. The sample blank value was also lowered if the serum volume was increased (Tab. 3). The significance of this effect, which could be caused by an inhibitor of the side reaction present in human serum, requires further investigation.

Further experiments were performed to investigate whether sera diluted with various solutions have higher activities than undiluted samples. Serum samples

Tab. 4

The activity of the glutamate dehydrogenase in undiluted (A) and serum sample diluted 1 + 2 with 0.9% NaCl (B), Qualtrol (C) and albumine solution (70 g/l, D). All values (U/l) are means from 2 determinations

Patient No.	A	B	C	D
1	26.8	29.4	28.2	31.8
2	26.8	28.2	26.1	28.2
3	30.0	34.2	31.8	35.7
4	16.4	17.4	17.4	19.2
5	15.5	13.8	12.0	16.5
6	25.7	24.0	21.9	27.3
7	22.8	27.3	23.4	30.6
8	8.0	11.1	8.7	9.6
9	24.7	24.0	23.1	27.8
10	13.5	13.8	11.2	14.7
mean	21.0	22.4	20.4	24.1

Tab. 5

The stability of the glutamate dehydrogenase activity in Enzatrol (Lot. No. ET 233 A) and human sera if stored at room temperature (A) and at 4°C (B)

Time of storage		2 h	5 h	24 h	48 h	5 days	10 days
Enzatrol* A		5.6 $\pm$ 0.2	5.2 $\pm$ 0.3	4.6 $\pm$ 0.4	4.2 $\pm$ 0.4	4.0**	3.4**
B		5.8 $\pm$ 0.2	5.5 $\pm$ 0.3	5.2 $\pm$ 0.3	4.8 $\pm$ 0.2	4.6**	4.1**
Pat. 1	B	7.2	7.2	7.2	7.3	—	—
Pat. 2	B	11.4	11.8	12.9	—	12.5	—
Pat. 3	A	19.9	19.2	19.2	19.1	19.0	—
B		19.9	19.9	19.9	19.8	19.8	—
Pat. 4	A	7.0	7.0	6.7	6.4	6.0	—
B		7.3	7.3	7.0	6.8	6.5	—
Pat. 5	A	30.8	—	28.2	28.2	28.2	17.0
B		30.8	—	29.5	29.5	29.4	23.7
Pat. 6	A	12.6	—	12.2	11.9	10.9	6.7
B		13.5	—	12.7	12.5	11.7	9.8

\* mean values from 7 (with standard deviation) and \*\* from 2 experiments

from 10 different patients were diluted 1 + 2 with 0.9% NaCl, Qualtrol and albumin solution. The glutamate dehydrogenase activity was determined in these preparations and in undiluted samples. The results compared in Table 4 indicate that Qualtrol is suitable for the dilution of serum samples with glutamate dehydrogenase activities above 30 U/l. Both Qualtrol and the albumin solution contained no detectable activity of glutamate dehydrogenase.

Serum should be stored no longer than 5 days at 4°C before the glutamate dehydrogenase activity is measured (Tab. 5). At the present time, Enzatzol is the only commercial available serum which is suitable for the quality control of the glutamate dehydrogenase determination. However, the activity of this enzyme declines constantly, especially during the first hours (Tab. 5). Therefore, we use Enzatzol for the control of precision 24 h after its reconstitution.

### References

1. SCHMIDT, E. & SCHMIDT, F. W. (1968), *Deut. Med. Wochenschr.* 93, 1153—1155. — 2. WAHLS, E. & ARNDT-HAUSER, A. (1968), *Ärzt. Lab.* 14, 397—407. — 3. SCHMIDT, E. & SCHMIDT, F. W. (1969), *Schweiz. Med. Wochenschr.* 99, 514—521. — 4. SCHMIDT, E. (1970), in: H. U. BERGMAYER *Methoden der enzymatischen Analyse*, Verlag Chemie, Weinheim 2, Aufl., 607—613. — 5. Standardisierung von Methoden zur Bestimmung von Enzymaktivitäten in biologischen Flüssigkeiten (1972), *this j.* 10, 182—192. — 6. HAECKEL, R. (1972), *this j.* 10, 235—242. — 7. KAISER, H. (1965), *Z. Analyt. Chemie* 209, 1—18. — 8. HAECKEL, R. & PORTH, A. J. (1972), *this j.* 10, 91—94. — 9. HAECKEL, R. & HAECKEL, H. (1972), in *Fortschritte der Klinischen Chemie, Enzyme und Hormone* (KAISER, E., ed.), Verlag der Wiener Medizinischen Akademie, Wien, 115—119.

Priv.-Doz. Dr. R. Haeckel  
3000 Hannover  
Karl-Wiechert-Allee 9