

J. Clin. Chem. Clin. Biochem.
Vol. 16, 1978, pp. 307–311

Determination of Perazine Serum Levels by Gas Liquid Chromatography under Clinical Routine Conditions

By J. Schley, E. Riedel, and B. Müller-Oerlinghausen

Institute of Biochemistry and Department of Psychiatry, Free University of Berlin

(Received September 13, 1977/January 18, 1978)

Summary: A quantitative gas liquid chromatographic method for the determination of serum levels of perazine (10-[3'-(1''-methyl-4''-piperazinyl)-propyl]-phenothiazine) is described. Perazine is used as a neuroleptic drug. The main problem consists in optimizing the chromatographic system. A sensitivity of appr. 60–150 nmol/l (20–50 µg/l) serum is achieved. Examples of optimization, analyses with patient samples, and the reproducibility of the results are presented.

Gasliquidchromatographische Bestimmung der Perazin-Konzentration im Serum unter klinischen Routinebedingungen

Zusammenfassung: Es wird eine quantitative gasliquidchromatographische Methode zur Bestimmung von Perazin (10-[3'-(1''-Methyl-4''-piperazinyl)-propyl]-phenothiazin) im Serum beschrieben. Perazin wird als Neuroleptikum eingesetzt. Das methodische Hauptproblem ist die Optimierung des chromatographischen Systems, wodurch eine Empfindlichkeit bis zu 59–148 nmol/l (20–50 µg/l) erreicht wird. Reproduzierbarkeit und praktische Anwendbarkeit unter klinischen Bedingungen werden dargestellt.

Introduction

The quantitative determination of perazine (10-[3'-(1''-methyl-4''-piperazinyl)-propyl]-phenothiazine¹⁾ in human serum samples under clinical routine conditions presents an analytical problem which has not yet found a satisfactory solution. Therapeutic serum levels of perazine can be expected to be within the range of approx. 59–885 nmol/l (20–300 µg/l) (1). Spectrophotometric and fluorometric methods are unsuitable, owing to the codetermination of metabolites (2, 3). The chromatographic separation of serum eluates is an indispensable prelude to the detection step, but this entails a loss of material and reduced sensitivity. The autoxidability of the phenothiazines provides additional analytical problems.

Breyer & Villumsen (1) have developed a direct densitometric thin layer chromatographic method, which has been shown to be very sensitive but also rather time-consuming (4). The use of gas liquid chromatography (GLC) entails some special difficulties, which are

described in detail below. The very sensitive electron capture detectors, which have been used for determination of phenothiazine derivatives such as chlorpromazine or perphenazine, cannot be applied in the case of perazine. For the determination of 20–50 ng perazine by use of the FID, however, the analytical procedure must be fully optimized. In the method published recently by Vanderheeren et al. (5) most of the steps are not optimized, so that determinations are limited to samples containing from 1.47 nmol (= 100 µg/l serum with 5 ml serum extracted) upwards. It is, however, possible, to perform 20 GLC determinations of perazine daily per column with a lower detection limit of 59 nmol/l (20 µg/l), providing the following optimization procedures are taken into consideration.

Material and Methods

Material

All reagents were Merck p.a. quality. Silicone XE 60, Silicone OV-12, dichlormethylsilane, Chromosorb G and Gaschrome Q 80/100; 100/120 mesh were ordered from WGA-Düsseldorf.

¹⁾ Taxilan, Promonta, D-2000 Hamburg 26

Procedure

Preparation of glass vessels

Cleaning of all glass vessels should be done without detergents. All glass vessels, after being cleaned with hot water, are boiled in about 6 mol/l HCl, rinsed with double distilled water, dried, and treated with 50 ml/l dichloromethane in toluene. They are dried finally at 120°C for 1 hour.

Extraction

5 ml serum (deep frozen at -15°C and thawed shortly before processing) with addition of 1 ml of 1 mol/l NaOH are kept in a boiling water bath for 15 min, then cooled down to room temperature and extracted 3 times with 5 ml *n*-heptane. The combined extracts are evaporated at room temperature or with light heating (not > 35°C) under nitrogen.

Analysis

GLC was performed with a Varian 1520B using flame ionisation detectors. Column: glass, 190 cm, 4 mm ϕ , silanized, packed with 3% silicone OV-17 on Gaschrom Q 100/120 mesh. Injection port: 270°C detector: 300°C. Oven: 270°C. Carrier gas flow 30 ml N₂/min. H₂: 30 ml/min. Air: 300 ml/min. Attenuation 2 or 4×10^{-11} . The evaporated heptane extracts are taken up with 25 μ l standard solution (internal standard: 48.46 nmol/l stigmasterol in ethanol); they must be kept cool at 0–4°C and tightly closed until injection. After careful stirring with a spatula 1 μ l is injected. The column is checked daily by running a calibration curve. The quantitative evaluation is performed by peak height comparison.

Results and Discussion

Extraction procedure

We investigated the distribution coefficients and elution characteristics of the following solvents: Methylene dichloride, chloroform, ethylene dichloride, benzene, toluene, *n*-heptane, *n*-heptane/isoamyl alcohol (985 ml + 15 ml) in relationship to the pH of the aqueous phase (4). *n*-heptane was found to be the most suitable solvent for the GLC procedure: at pH 10.4 the distribution coefficient was 5.7 and the coextraction of serum factors was the lowest with regard of the temperature range used for perazine detection (fig. 1 a). Sufficiently low solubility of perazine base in water is achieved in an alkaline milieu, i.e. pH 8.5. Short heating of the alkalized serum samples in a boiling water bath results in a better separation of the organic and the aqueous phase without loss of the drug (tab. 1). Inactivation of the surface of the glass vessels e.g. by silanization leads to considerably increased recovery values. The use of detergents for cleaning of the test tubes etc. must be avoided under any circumstances. Precleaning is done best by boiling with about 6 mol/l hydrochloric acid.

Optimization of the GLC system

As stationary phases we have examined silicone-XE-60 and silicone-OV-17 (3%) on Chromosorb G or Gaschrom Q (80/100 or 100/120 mesh resp.). Resolution with the 80/100 mesh supports was lower and the loss of perazine greater than with the 100/120 mesh ab-




Fig. 1. Gas chromatogram (different series at different days).

1a: Blank serum

1b: Perazine/stigmasterol in ethanolic solution on 3% silicone-OV-17 (Gaschrom Q, 100/120 mesh)
Attenuation 8×10^{-12} As

1c: Serum of a patient treated with perazine
Attenuation 8×10^{-12} As

Tab. 1. Recovery of perazine extracted from aqueous standard solutions or from spiked human sera under different analytical conditions.

Perazine concentration [nmol/l]	Recovery rate			
	From aqueous solution		From spiked sera	
	No treatment	Heated 15 min at 100°C	Heated 15 min at 100°C	Unsilanized glassware
	N [%]	N [%]	N [%]	Silanized glassware
147.5	10 94	8 90	5 52	8 76
295	15 86	7 89	6 53	6 73
590	16 88	6 89	4 59	6 72

sorbents. The most suitable stationary phase proved to be 3% silicone-OV-17 on Gaschrom Q 100/120 mesh. Optimization is achieved by determination of the *Van Deemter* hyperbolic function which relates the column efficiency – expressed as height equivalent to a theoretical plate (HETP) – to carrier gas flow changes; in addition, peak height, peak area and retention time were taken into consideration (fig. 2). Measurements were performed with standard solutions of 2.95 mmol/l perazine in ethanol. The aim is to find the carrier gas flow rate at minimal retention time (in




Fig. 2. Column efficiency (height equivalent to a theoretical plate HETP \circ — \circ), peak height (\circ — \circ), peak area (\circ — \circ) and retention time (\times — \times) as a function of carrier gas flow (N_2). Analyses were performed with 2.95 nmol perazine. Technical data ref. text.

order to increase the laboratory capacity) and at maximum peak height or area (in order to increase the sensitivity). From figure 2 it can be seen that it is not reasonable to make use of the highest possible separation performance (achieved at $u \leq 10$ ml/min), because then the peak height is too small and retention time too long. On the other hand, it does not make much sense to increase the carrier gas flow beyond $u = 30$ ml/min, because this does not result in a substantial reduction of the retention time, but in a considerable loss of sensitivity. Figure 2 shows that the optimal carrier gas flow for highest sensitivity and shortest retention time is 30 ml/min. Another variable to be optimized is the column temperature. With increasing column temperature, peak heights are elevated and retention times shortened (fig. 3), thereby leading to an improved sensitivity and increased analytical capacity. If the full sensitivity of the chromatographic system is used the perazine

peak comes within the tailing of the solvent peak (at column temp. $> 250^\circ\text{C}$); electronic integration techniques are, therefore, not recommended. The quantitative evaluation should be done by comparison of the peak heights with those of an internal standard.

Reproducibility

Once the optimal values for the carrier gas flow and the column temperature have been determined, their daily control is essential for routine analytical work. This can be best achieved by setting up a calibration curve before every course of analysis. 29, 59, and 118 nmol (= 10, 20, or 40 μg , resp.) perazine were used with 48.46 nmol/l stigmasterol as internal standard (fig. 4). The evaluation by comparison of the peak heights does not give a linear calibration curve in all cases. The correct ratio between perazine and stigmasterol concentration in the solvent must be carefully checked. Evaporation of the solvent leads to a higher concentration and thereby to changes of the peak height ratio which may result into evaluation errors. Figure 5 shows the decrease of the peak height ratios with increasing concentration of perazine and stigmasterol in the ethanolic solution (at always the same concentration ratio 1:1). Therefore, extracts which have been taken up with the standard solution should be kept cool (0 – 4°C) and analysed the same day.

The recovery and the analytical error of analyses obtained with perazine-spiked blank sera can be seen from table 3. The quantitative evaluation of perazine determinations in serum is done by referring to the daily calibration curve and multiplication with the correction factor 1.35. The reproducibility of the method was also




Fig. 3. Retention time (—) and peak height (—) of perazine (\bullet) and stigmasterol (\circ) as a function of the column temperature. Technical data ref. text.

Tab. 2. Optimization of the GLC system.

Technical data ref. text. Analyses were performed with 2.95 nmol perazine.

u = carrier gas flow; t_R = retention time

P = peak height; A = peak area; HETP = $1000 l/N$

(l = column length (m); N = number of theoretical plates

= $16 \left(\frac{t_R}{W} \right)^2$ with W = peak width (min).

All values are means from 3–10 single measurements.

u [ml/min]	t_R [min]	P [mm]	A [cm ²]	HETP [mm]	N
10	10.0	21	3.4	0.38	5000
22	7.0	100	5.6	0.48	4000
30	5.8	101	5.4	0.54	3500
36	5.0	96	4.5	0.60	3200
40	4.3	50	2.2	0.62	3100
47	4.3	40	1.6	0.56	3400




Fig. 4. Perazine calibration curve using peak height comparison with stigmasterol as internal standard.




Fig. 5. Peak height ratio (ref. fig. 4) as a function of the concentration of perazine and stigmasterol (1:1) in the ethanolic solution.

Tab. 3. Recovery of perazine after extraction from spiked human blank sera.

Perazine [nmol/l]	N	Recovery [nmol/l] $\bar{x} \pm s_x$	Re- cov- ery frac- tion	Correc- tion factor
147.50	50	8	112.10 ± 14.36	0.76
295.00	100	17	216.91 ± 36.81	0.74
590.00	200	6	419.02 ± 51.62	0.72

examined by analysing the same serum samples twice, i.e. 14 samples which were split immediately after collecting the blood and processed independently of each other within an interval of 5–8 days. The mean difference

Tab. 4. Reproducibility of perazine determinations in serum samples of 14 out-patients. Each sample was analyzed twice (I–II) within 5–8 days.

Patient	Perazine I [nmol/l]	Perazine II [nmol/l]	Patient	Perazine I [nmol/l]	Perazine II [nmol/l]
D. M.	100.30	103.25	W. L.	218.30	183.90
K. J.	171.10	162.25	B. H.	76.70	73.75
B. G.	227.15	212.40	R. J.	516.25	501.50
K. B.	182.90	171.10	S. F.	289.10	265.50
N. D.	398.25	377.60	S. H.	197.65	182.90
K. W.	454.30	442.50	Z. A.	286.15	271.40
H. U.	702.10	663.75	P. H.	59.00	44.25

Tab. 5. Retention time (t_R), peak height (P) and relative peak height "P" as compared to perazine.

Compound	t_R [min]	P [mm]	"P"
Perazine	6.0	198	1.00
Desmethylperazine	7.5	75	0.37
Perazine sulfoxide	23.0	30	0.15

between the double determination was 17.38 nmol/l corresponding to a ratio of 0.065 between mean difference and mean value (tab. 4). The precision within one day was determined as $VK = 2.55\%$ ($N = 6$; $\bar{x} = 364.3$ nmol/l; different extractions of one pooled sample). Relative retention time was 0.59.

Specificity

The retention time of the following psychotropic compounds were determined and it was found that they do not interfere with the perazine measurement: butaperazine, trifluoperazine, levomepromazine, chlorpromazine, haloperidol, imipramine, amitriptyline, nortriptyline.

Concluding remarks

The lowest detection limit of the method as described above is appr. 60 nmol/l. A clinical study which should demonstrate the practicability of the method under routine conditions, i.e. in perazine treated out-patients, has just been finished (6). Two main metabolites, desmethylperazine and perazine sulfoxide cannot be determined together with perazine in human sera. Perazine, desmethylperazine, and perazine sulfoxide are well separated from each other by this GLC procedure; however, the retention time of desmethylperazine is the same as that of a serum lipid peak, and the sensitivity of the system for the sulfoxide is about 6 times lower than for perazine.

An independent GLC method for the determination of perazine in body fluids is thus available, in addition to the TLC method quoted above, thus providing the prerequisite for a systematic comparative control of these two clinical methods. Such internal comparative control by two independent methods appears to us much more relevant than the quality control by multicentered studies, which may result in a better comparability of the results of different laboratories, but which cannot exclude systematic analytical errors in the case of methods not yet standardized. The comparative study of Richens (7) provides much evidence for this argument. As the determination of neuroleptics or antidepressant drugs in human sera requires at least a 100 times higher

analytical sensitivity than in the case of anticonvulsant drugs, a rigid and comparative control of the methods is essential.

Acknowledgements

Thanks are due to Dr. U. Breyer for providing us with samples of desmethylperazine and perazine sulfoxide. Perazine base was a gift from Promonta, Hamburg.

This work was supported by a grant from the Deutsche Forschungsgemeinschaft („Klinische Pharmakologie“; He 916/1-2).

References

1. Breyer, U. & Villumsen, K. (1976), *Eur. J. Clin. Pharmacol.* **9**, 457-465.
2. Hadass, H., Kluwe, S. & Fähndrich, Ch. (1972), *Pharmakopsychiatr.* **5**, 53-69.
3. Breyer, U., Müller-Oerlinghausen, B. & Mauruschat, W. (1976), in *Psychotherapeutic Drugs* (Forrest, I. S. & Usdin, E., eds.), Marcel-Dekker, New York.
4. Müller-Oerlinghausen, B., Riedel, E. & Schley, J. (1977), *Int. J. Clin. Pharmacol.* **15**, 366-372.
5. Vanderheeren, F. A. J., Theunis, D. J. C. J. & Rosseel, M. T. (1976), *J. Chromatogr.* **120**, 123-128.
6. Pietzcker, A., Schley, J. & Müller-Oerlinghausen, B. (1977), *Arzneim.-Forsch. (Drug Res.)*, in press.
7. Richens, A. (1975), in *Clinical Pharmacology of Anti-Epileptic Drugs* (Schneider, H. et al., eds.), pp. 293-302, Springer-Verlag, Berlin.

Prof. Dr. E. Riedel
Institut für Biochemie
der Freien Universität Berlin
Limonenstr. 7
D-1000 Berlin 45

