

J. Clin. Chem. Clin. Biochem.
Vol. 18, 1980, pp. 209–214

Quantitative Determination of Clobazam in Serum and Urine by Gas Chromatography, Thin Layer Chromatography and Fluorometry

By P. Hajdú, M. Uihlein and D. Damm

Hoechst Aktiengesellschaft, 6230 Frankfurt/Main 80

(Received August 24/October 18, 1979)

Herrn Prof. Dr. Rolf Sammet zu seinem 60. Geburtstag gewidmet


Summary: The procedures available for determination of clobazam (Frisium, Hoechst) are gas chromatography, fluorometry, and thin-layer chromatography. The study presents detailed descriptions of analytical procedures appropriate for routine determinations in serum and urine, and results from human trials. Moreover, the physico-chemical properties of clobazam, viz., solubility, distribution, and protein binding are given.

Bestimmung von Clobazam in Serum und Urin mittels Gaschromatographie, Dünnschichtchromatographie und Fluorometrie

Zusammenfassung: Die vorliegende Arbeit gibt eine ausführliche Beschreibung von gaschromatographischen, dünn-schichtchromatographischen und fluorometrischen Analysemethoden zur routinemäßigen Bestimmung von Clobazam (Frisium, Hoechst) und seinen wichtigsten Metaboliten in Serum und Urin. Es werden Ergebnisse von Untersuchungen am Menschen präsentiert und darüber hinaus Angaben zu folgenden physikalisch-chemischen Eigenschaften des Präparats gemacht: Löslichkeit, Verteilung und Protein-Bindung.

Introduction

Methods are described for the determination of the new anxiolytic agent clobazam¹⁾ (a 1,5-benzodiazepine (1)) in serum (gas chromatography and fluorometry) and urine (thin-layer-chromatography). Clinically relevant physico-chemical data are also presented.



Gaschromatographic Determination in Serum

Clobazam and its main metabolite in serum, N-des-methylclobazam, are determined by gas chromatography using procedure already described for 1,4-benzodiazepines (2–4).

¹⁾ Frisium®, Hoechst Aktiengesellschaft.

Reagents

Ethyl ether, ethyl acetate, both reagent grade; internal standard: 1 mg/l diazepam in demineralized water.

Equipment

Glass stoppered centrifuge tubes volume about 10 ml; glass stoppered test tubes with tapered end after *Beckett*; Vortex mixer; refrigerating centrifuge; gas chromatograph; integrator HP 3380.

Processing

To 1 ml serum in a centrifuge tube add 0.1 ml of the internal standard solution, and mix with 5 ml ethyl ether on a Vortex mixer for 30 s. Centrifuge for 5 min in a refrigerating centrifuge at 0–5°C, transfer about 4 ml of the organic phase to a *Beckett* tube, and evaporate to dryness at about 40 °C under nitrogen flow. Take up the residue with 50 µl ethyl acetate, and keep in a refrigerator until analysis.

Gas-chromatography

Apparatus:	Hewlett-Packard 5700 with nitrogen-specific FID 18 789 A and glass-jet	
Column:	Glass column 1.20 m × 4 mm	
Stationary phase:	2% OV 101 on Chromosorb W/AW-DMCS, 120–150 µm	
Temperature:	Oven:	240 °C isothermal
	Injection port:	300 °C
	Detector:	300 °C

Carrier gas: Helium
 Flow rate: 40 ml/min
 Sample volume applied: 1–2 μ l
 Duration of analysis: about 6 min

Figure 1 shows the chromatograms of a serum blank, a serum blank with the addition of 0.05 mg/l of clobazam and N-desmethyclobazam, and a patient's serum containing 0.2 mg/l clobazam and 0.11 mg/l N-desmethyclobazam. Evaluation was carried out with the integrator by an internal standard procedure, in which the peak areas of the compounds to be determined are compared with the peak area of diazepam.

Evaluation may also be performed by comparing the corresponding peak heights in the chromatogram. The linearity of this evaluation method was verified with mixtures of the pure compounds in ethyl acetate in the range 1 to 25 ng injected ($r = 0.9992$).

Precision and accuracy

Mixtures of clobazam and N-desmethyclobazam were diluted in serum. The serum samples were split in order to obtain six identical dilution series.

These series were analysed consecutively and gave the following results (tab. 1).

Tab. 1. Gas-chromatographic determination of clobazam and N-desmethyclobazam in serum.

Added (mg/l)	Found Clobazam (mg/l)		N-Desmethyclobazam mg/l	
0.5	0.51 ± 0.012	(2.4%)	0.53 ± 0.016	(3.0%)
0.2	0.21 ± 0.008	(3.8%)	0.21 ± 0.028	(13.0%)
0.1	0.100 ± 0.0033	(3.3%)	0.091 ± 0.009	(9.9%)
0.05	0.055 ± 0.0026	(4.7%)	0.052 ± 0.009	(17.3%)
0.02	0.019 ± 0.0008	(4.2%)	0.020 ± 0.002	(10.0%)
0.01	0.011 ± 0.0010	(9.1%)	—	—
Blank*	0.007 ± 0.0017		0	
Bias	$+ 0.004 \pm 0.005$ mg/l		$+ 0.007 \pm 0.015$ mg/l	

*The blank value was obtained from a compound interfering with clobazam (retention time = 3.59 min versus 3.52 min for clobazam). As the separation of these compounds is extremely difficult and the blank value is only in the order of magnitude of the detection limit, we suggest a correction of the results.

For clobazam, the relative error of the method increases with decreasing values, which is usually the case in chromatographic trace determinations. The absolute errors correlate well and form a straight line, which does not go through the origin. This permits us to describe the error of the method as the sum of a constant and relative error. Thus, the precision of the method calculated for clobazam amounts to 0.002 mg/l + 2.2% of the value measured.




Fig. 1. Determination of clobazam and N-desmethyclobazam in serum by gas chromatography.

- a) Serum blank
- b) Admixture of 0.05 mg/l each to serum blank
- c) Serum levels after 5th administration of clobazam

Glass column 1.2 m x 4 mm 2% OV 101 on chromosorb W/AW-DMCS 240 °C isothermal.

I: Internal standard

C: Clobazam

N: N-desmethyclobazam

U: Unknown

Although the above-mentioned correlation is not satisfactory for N-desmethyclobazam, we prefer to describe the error dependency in the same way, resulting in $0.006 \text{ mg/l} + 2.5\%$ of the value measured.

The precision of the method thus obtained enables the limit of detection for quantitative determinations to be easily assessed. If the limit of detection is defined as $t_p \times \text{S. D.}$ (the concentration being 0), the following values are obtained: ($n = 6$, $t_p: 0.05 = 2.06$)

0.005 mg/l for clobazam

0.015 mg/l for N-desmethyclobazam

Selectivity

The method measures clobazam and N-desmethyclobazam selectively. The determination of N-desmethyclobazam in therapeutical levels is only possible if a glass column and FID with glass-jet are used. The following synthesized polar metabolites can certainly not be detected in therapeutical levels:

4'-hydroxyclobazam

4'-(3')-hydroxy-3'(4')-methoxyclobazam

4'-hydroxy-N-desmethyclobazam

4'-(3')-hydroxy-3'(4')-methoxy-N-desmethyclobazam

Results obtained from a multiple-dose human study

In a human study, clobazam was administered to seven volunteer test persons for 10 days daily at 9 a.m. in a dose of 10 mg and at 7 p.m. in a dose of 20 mg. The blood samples were taken repeatedly on the first day, and before each administration on the following days. (5).

On the first day, the maximum serum level of clobazam (0.31 mg/l) was reached two hours after the administration. The elimination half-life was eight hours. N-desmethyclobazam was detected in the blood not earlier than eight hours after administration in a concentration of only 0.05 mg/l.

On the fourth day of treatment, the clobazam levels reached a steady state of $0.90 \pm 0.23 \text{ mg/l}$. N-desmethyclobazam which had a half-life of 50 hours reached a steady state approximately on the 9th day of treatment at serum levels of $3.6 \pm 1.3 \text{ mg/l}$. An example for the course of serum levels is shown in figure 2.

Fluorometric Determination in Serum

Clobazam can also be determined by fluorometry. It does not show native fluorescence, but can be transformed into a fluorophore by irradiation with short-wave UV light (6). The advantage of the method is that it is independent of a blank value as the sample is measured before and after irradiation. This transformation shows good reproducibility (cf. 1. c. (7)).




Fig. 2. Concentrations in human serum following multiple dose of clobazam (30 mg/day).

● clobazam
○ N-desmethyclobazam

Reagents

Ethyl ether, ethyl alcohol, both reagent grade; Standard solution: 0.5 mg/l clobazam in ethyl alcohol.

Equipment

Glass stoppered centrifuge tubes, volume about 10 ml; Vortex mixer; centrifuge; Pen-Ray lamp 11 Sc 1 (Ultra-Violet Products, Inc., San Gabriel, California); irradiation chamber with two 1 cm cells; fluorometer Perkin Elmer MPF 3.

Processing

To 1 ml serum in a centrifuge tube add 5 ml ethyl ether, and shake for 30 s on a Vortex mixer. Centrifuge briefly, pipet off 4 ml from the supernatant ether phase into a test tube, and evaporate to dryness at about 40°C under nitrogen flow. Take up the residue with 2 ml ethyl alcohol (sample solution) and keep protected from light until fluorometric measurement.

Fluorometric measurement

Adjustment:

Excitation wavelength	364 nm
Excitation slit	3 nm
Emission wavelength	400 nm
Emission slit	6 nm
Sample sensitivity	30

Calibration:

In the irradiation chamber, irradiate 2 ml standard solution with UV light in a 1 cm cell for about one minute. The irradiation time depends on the intensity of the lamp and has to be readjusted from time to time. Adjust the scale reading obtained from the irradiated standard solution to 100% amplitude by means of the sample adjust button. This adjustment remains unchanged during sample measurement.

Fill the sample solution into a cell and measure (value 1). Subsequently, irradiate the sample solution for exactly the same time as the standard solution and measure again (value 2).

Value 1 (non-irradiated solution) is regarded as blank value and subtracted from value 2 prior to calculation. Calibration curves are established with sera containing standard quantities of the compound. Calibration is linear up to a concentration of 5 mg/l.

Precision and limit of detection

The data (tab. 2) show a constant precision of the method of 0.016 mg/l within the concentration range 0.1 to 1.0 mg/l serum. Consequently, the limit of detection (cf. paragraph on gas chromatography) is 0.04 mg/l.

Tab. 2. Results obtained from admixtures to serum.

Added (mg/l)	Found (mg/l)	n
0.1	0.098 ± 0.013	5
0.2	0.200 ± 0.017	15
0.3	0.310 ± 0.020	5
0.4	0.405 ± 0.013	15
0.5	0.500 ± 0.016	5
0.6	0.596 ± 0.023	10
0.8	0.794 ± 0.018	10
1.0	1.005 ± 0.010	10
Blank value	0 ± 0.006	17
Bias	+ 0.001 ± 0.005 mg/l	

Selectivity

We examined whether the metabolites mentioned above show fluorescence under the same conditions as clobazam: the fluorescence of N-desmethyloclobazam amounts only to one third of that of clobazam. 4'-Hydroxyclobazam and 4'-hydroxy-N-desmethyloclobazam show one tenth of the fluorescence of clobazam. The two hydroxy-methoxy metabolites show no fluorescence. 1,4-Benzodiazepines do not interfere with the method.

Although the fluorometric method is less specific than gas chromatography, it is a useful alternative for single-dose serum kinetics (cf. next paragraph), where these metabolites are only found in negligible serum concentrations (8). Stewart et al. showed (7) that the selectivity can be improved by using hexane as extraction agent.

Comparison of the methods

The quality of both methods is best demonstrated by parallel determinations in human trials. The results obtained from parallel determinations performed 0 to 24 hours after administration to 10 test persons correlate well (cf. fig.3);

$$n = 60$$

$$r = 0.967$$

$$C_{GC} = 0.96 \times C_{FI} + 0.015 \text{ mg/l}$$

$$S_{y,x} = S_{x,y} = 0.042 \text{ mg/l}$$

It should be mentioned that this correlation exclusively refers to clobazam, because N-desmethyloclobazam was found only in negligible amounts.

Thin-layer chromatographic determination in urine

Clobazam and the main metabolites in man, dog, and monkey – N-desmethyloclobazam, 4'-hydroxyclobazam, and 4'-hydroxy-N-desmethyloclobazam – can also be quantitatively determined by thin-layer chromatography. The compounds are easily separated on silica gel, using the mobile phase chloroform/n-heptane/methanol (85 ml + 10 ml + 5 ml). Quantitative evaluation is carried




Fig. 3. Correlation of concentrations of clobazam in human serum obtained by gas chromatography and fluorometry.

out by reflectance measurement with the usual TLC evaluating devices (densitometers) at 230 nm.

As a large amount of the metabolites is excreted as conjugates (8), the urine has to be treated with glucuronidase before determination.

Reagents

Ethyl acetate, chloroform, *n*-heptane, methanol, all reagent grade; acetate buffer pH 5.5 (24 g CH₃COONa · 3 H₂O + 240 ml 0.1 mol/l acetic acid/l H₂O); phosphate buffer pH 7.5 (15.3 g Na₂HPO₄ · 2H₂O + 1.9 g KH₂PO₄/l H₂O saturated with NaCl); 0.01 mol/l sodium hydroxide solution; 0.01 mol/l hydrochloric acid; β-glucuronidase (Boehringer, Mannheim); ready-prepared TLC plates Si 60 F₂₅₄ 20 × 20 cm (silica gel with fluorescence indicator, E. Merck, Darmstadt)

Equipment

Glass stoppered centrifuge tubes, volume about 10 ml; Vortex mixer; water bath; centrifuge; Desaga mobile phase trough; Desaga autospotter; TLC-spectrophotometer PMQII with micro-optics (Carl Zeiss); recorder Servogor 310.

Processing

a) Free compounds

To 1 ml urine in a centrifuge tube add 1 ml phosphate buffer pH 7.5, and mix with 5 ml ethyl acetate on a Vortex mixer for 30 s. Centrifuge for 5 min at 4,000 rpm, transfer about 4.5 ml of the organic phase to a centrifuge tube, treat with 1 ml 0.01 mol/l NaOH, and shake again on a Vortex mixer for 30 s. Centrifuge, discard the (lower) aqueous layer, and treat about 4.2 ml of the organic phase with 1 ml 0.01 mol/l HCl in the same way. Subsequently, evaporate 4 ml of the organic phase at about 40°C under nitrogen flow. Take up the residue with 100 μl methanol before application.

b) Conjugated compounds

Treat 1 ml urine in a centrifuge tube with 1 ml acetate buffer pH 5.5 and 5 ml ethyl acetate on a Vortex mixer for 30 s in order to separate the free compounds. Extract again with 5 ml ethyl acetate. Treat the remaining aqueous phase²⁾ with 10 μl

²⁾ Addition of calibration compounds is made at this stage, since the appropriate glucuronides are not available.

β -glucuronidase, and digest at 37°C in a water bath for 48 h. Subsequently add 1 ml phosphate buffer pH 7.5 and proceed as described under "Free compounds".

Application and development

The autospotter is able to apply up to 20 samples simultaneously onto a TLC plate, whereby optimal adjustment secures minimal spot size.

Apply onto a plate 14 extracts of samples to be measured, 5 extracts of known admixtures, and a sample of pure substance.

Adjustment:

Application speed	3
Heating	3
Ventilation	3
Application volume	50 μ l

The plates are developed simultaneously in a mobile phase of chloroform/*n*-heptane/methanol (85 ml + 15 ml + 5 ml) at room temperature.

Length of run: 16 cm

Development time: about 100 min

R_f values:

Clobazam	0.90
N-Desmethyclobazam	0.59
4'-Hydroxyclobazam	0.48
4'-Hydroxy-N-desmethyclobazam	0.31

Measurement

Adjustment of TLC spectrophotometer PMQ II:

Slit length 2 mm, slit width 14 mm (equivalent to 0.5 \times 3.5 mm with microoptics)

Wavelength 230 nm

Speed of plate reek 30 mm/min

Adjust the plate reek to the first substance spot and scan the plate within the corresponding R_f range.

Evaluation

Evaluation on the basis of the peak areas is done by hand or by means of an integrator. For each plate, the calibration is calculated separately from admixtures to urine. Linearity is secured up to a concentration of 1 mg/l urine. At higher concentrations, either less urine is used or a smaller aliquot is applied onto the plate.

Results from addition of standard substances to urine

Since the original substance is not found in urine (6), only the metabolites were added. These dilutions were split into 21 identical samples and processed as described for free and conjugated compounds. The following results were obtained (tab. 3).

Free compounds

Precision for all three metabolites is 0.025 mg/l + 4% of value. The limit of detection calculated from this is 0.05 mg/l (c. f. paragraph on gas chromatography).

Conjugated compounds

	Precision	Detection limit
N-desmethylclobazam	0.035 mg/l + 3% of value	0.07 mg/l
4'-Hydroxyclobazam	0.065 mg/l + 1.4% of value	0.13 mg/l
4'-Hydroxy-N-desmethylclobazam	0.025 mg/l + 2% of value	0.05 mg/l

Tab. 3. Thin-layer chromatographic determination of metabolites of clobazam in urine.

Added (mg/l)	Found N-Desmethylclobazam (mg/l)	Found 4'-Hydroxyclobazam (mg/l)	Found 4'-Hydroxy-N-desmethyclobazam (mg/l)
a) Processed as free compounds			
1.0	0.99 \pm 0.06	1.00 \pm 0.06	0.99 \pm 0.06
0.5	0.48 \pm 0.05	0.51 \pm 0.05	0.53 \pm 0.06
0.2	0.20 \pm 0.05	0.22 \pm 0.03	0.24 \pm 0.03
0.1	0.09 \pm 0.01	0.10 \pm 0.02	0.09 \pm 0.03
Bias	- 0.01 \pm 0.01	+ 0.01 \pm 0.01	+ 0.01 \pm 0.03
b) Processed as conjugated compounds			
2.0	2.00 \pm 0.06	1.96 \pm 0.09	1.96 \pm 0.09
1.0	1.01 \pm 0.06	1.08 \pm 0.09	1.03 \pm 0.07
0.5	0.53 \pm 0.04	0.56 \pm 0.08	0.53 \pm 0.07
0.2	0.18 \pm 0.03	0.28 \pm 0.07	0.20 \pm 0.03
0.1	0.09 \pm 0.02	0.12 \pm 0.06	0.09 \pm 0.03
Bias	0 \pm 0.02	+ 0.04 \pm 0.05	0 \pm 0.03

Figure 4 shows the cumulative urinary excretion of the conjugated metabolites after a single oral dose of 40 mg clobazam given to a healthy volunteer.




Fig. 4. Cumulative excretion in urine after a single dose of 40 mg clobazam.

Subj. no. 1			
N-desmethyclobazam			
○—○	A0 = 0.032	K1 = 2.130	
	K2 = 0.360	T0 = 0.019	
4'-hydroxyclobazam			
x—x	A0 = 0.017	K1 = 4.063	
	K2 = 0.099	T0 = 0.001	
4'-hydroxy-N-desmethyclobazam			
●—●	A0 = 0.038	K1 = 1.767	
	K2 = 0.206	T0 = 0.001	

Physico-Chemical Data of Clobazam

Besides the determination methods, the following properties are of clinical interest:

UV absorption spectrum

The UV spectrum (fig. 5) was taken with a recording spectrophotometer (Beckman Acta M VI) at a concentration of 10 mg/l in ethanol. It shows a maximum at 228 nm ($A_{1\text{ cm}}^{1\text{ g/dl}} = 1450$) and a weak double maximum at 285 and 295 nm ($A_{1\text{ cm}}^{1\text{ g/dl}} = 75$).

Solubility

The solubility in water is 80 mg/l within the pH range 3 to 11.




Fig. 5. Absorption spectrum of clobazam (10 mg/l in ethanol) 1 cm optical pathway.

Distribution

The compound is freely soluble in lipids. Aqueous buffers containing 5 mg/l of clobazam and N-desmethyloclobazam (5 mg/l) were extracted with octanol. Within the pH range 1.2 to 12.4, both compounds were found practically quantitatively in the organic phase at the distribution equilibrium.

Protein binding

Protein binding of the compound was determined by equilibrium dialysis according to Scholtan (9). Solutions of 0.125–2.0 mg/l clobazam in human serum were dialyzed against 0.067 mol/l phosphate buffer pH 7 at room temperature overnight. At equilibrium clobazam was determined fluorometrically in serum and the corresponding buffer solutions. The adsorption isotherm after Freundlich (cf. fig. 6) gave the following equation:

$$\log c_{\text{bound}} = 0.925 \times \log c_{\text{free}} + 0.780 \pm 0.045$$

i.e., 85 to 90 percent clobazam is bound to serum protein at therapeutical levels (0.1–1.6 mg/l). Similar results were obtained with ultracentrifugation (Centriflo membrane cones from Amicon, Lexington, Mass./USA).




Fig. 6. Protein-binding of clobazam in human serum.
 $r = 0.996$

References

1. Barzaghi, F., Fournex R. & Mantegazza, P. (1973) *Arzneim.-Forsch.* 23, 683–686.
2. Marcucci, F. (1968) *J. Chromatogr.* 37, 318–320.
3. de Silva, J. A. F. & Puglisi, C. V. (1970) *Anal. Chem.* 42, 1725–1736.
4. Zingales, I. A. (1971), *J. Chromatogr.* 61, 237–252.
5. Rupp, W., Badian, M., Christ, O., Hajdú, P., Kulkarni, R. D., Taeuber, K., Uihlein, M., Bender, R. & Vanderbeke, O. (1979) *Br. J. Clin. Pharmacol.* 7, 51–57.
6. Hajdú, P. & Damm, D. (1976), *Arzneim.-Forsch.* 26, 2141–2145.
7. Stewart, J. T., Honigberg, I. L., Tsai, A. Y. & Hajdú, P. (1979) *J. Pharmaceut. Sci.* 68, 494–496.
8. Volz, M., Christ, O., Kellner, H.-M., Kuch, H., Fehlhaber, H.-W., Gantz, D., Hajdú, P. & Cavagna, F. (1979) *Br. J. Clin. Pharmacol.* 7, 41–50.
9. Scholtan, W. (1964) *Arzneim.-Forsch.* 14, 146–149.

Dr. P. Hajdú
Hoechst Aktiengesellschaft
Postfach 80 03 20
6230 Frankfurt 80